

Disc Grinder HK200

Mobile but sturdy sample preparation machine for spectrometer samples

- Mobile desk-top model
- For use with self-adhesive grinding paper
- Solid steel cabinet, cast iron disc pan and motor base
- Easy to clean
- Disc pan with exhaust outlet.

Technical specifications:

Power supply	230 V / 50 Hz, 370 W
Instrument dimensions (W x D x H)	270 x 290 x 280 mm
Alternative power	230 V / 60 Hz or 115 V / 60 Hz
Mass	22 kg
Revolution speed	2880 rpm
Grinding disc diameter	205 mm

Consumables and accessories

Self adhesive grinding papers

- 200 mm available as Al_2O_3 , SiC, ZrO_2 in various grades.

Sample holders

- Magnetic (square)
- Magnetic (round)
- Mechanical (round).

